

- Perfect Illumination
- Easy to Upgrade
- Affordable



Superior Light

Go LED in seconds

- Simply exchange the halogen bulb in your KaVo MULTIflex LUX[™] coupler with the KaVo MULTI LED
- No downtime or repair time: You can independently exchange the LED in your practice
- You have the flexibility to buy the KaVo MULTI LED or the KaVo MULTIflex™ LED couplers to suit your unique need.



Easily Retrofits

Halo-free illumination

- Color temperature similar to daylight (5500 deg K)
- Optimally focused 25,000 lux at bur-tip
- Uniform pattern prevents glare from scattered light



Lower Cost of Ownership

2

Lasts longer than traditional bulbs

- You do not have to buy new LED MULTIflex couplers
- You only need to purchase one KaVo MULTI LED per coupler
- The service life of the KaVo LED is several times that of the quartz halogen bulb



MULTI LED™ Bulb

KaVo MULTI LED for your MULTIflex couplers. KaVo's innovative solution is an LED that can easily replace your existing halogen lamp. MULTI LED offers efficient and superior light, with color temperature similar to daylight (5,500K and 25,000 lux) at the bur tip.

MULTIflex[™] 460 LED

6-pin fiber optic coupler with MULTI LED bulb in coupler. Anti-retraction valve. (Same as 465LRN but without water spray regulator.) Use with LCM Illuminator and tubing. Not sterilizable.



MULTIflex[™] 465 LED

6-pin fiber optic coupler with MULTI LED bulb in coupler. Spray regulator. Antiretraction valve. Use with LCM Illuminator and tubing. Not sterilizable.



Product	Description	KaVo Part No.
KaVo MULTI LED	MULTI LED bulb with daylight color temperature and perfect illumination	1007.5372
MULTIflex 460 LED	Six-pin coupler with KaVo MULTI LED	1007.3201
MULTIflex 465 LED	Six-pin coupler with KaVo MULTI LED bulb and spray regulator	1007.3197
MULTIflex 460 LE	Six-pin coupler with halogen bulb (upgradeable to KaVo MULTI LED)	1001.7599
MULTIFlex 465 LRN	Six-pin coupler with halogen bulb (upgradeable to KaVo MULTI LED) and spray regulator	0553.1550

KaVo MULTIflex coupling for every application

KaVo's innovative solution is an LED that can easily replace your existing halogen lamp. MULTI LED offers efficient and superior light, with color temperature similar to daylight (5,500K and 25,000 lux) at the bur tip. It is an economical way to upgrade your existing MULTIflex couplers to LED.



KaVo. Dental Excellence. 888-ASK-KAVO · www.kavousa.com

* Non LUX Instruments